

Optimized cell culture

Establishing protocols for automated harvesting of various cell lines from multiwell plates

Introduction

Culturing of various stable cell lines, including bacterial, yeast, hybridoma, insect and mammalian lines, and also stem cells, is an important task in the pharmaceutical industry. The cells need to be maintained and expanded for use in cell-based assays, such as siRNA screening, ADME and high content screening, as well as for clone selection. In addition to cell culture flasks (T-Flask, RoboFlask™), cells may be grown in multiwell plates – such as 6-well, 12-well and 24-well plates. This application note describes an automated procedure for harvesting commonly used human and animal cell lines that have been cultivated in multiwell plates, relieving laboratory personnel from tedious cell harvesting work and providing more time for other valuable tasks.

Material and Methods

Robotic Equipment

Experiments were conducted under sterile conditions; the workstation was housed in a biological safety cabinet and the liquid system thoroughly sterilized. Automated harvesting was performed using a Freedom EVO® workstation (Figure 1), equipped with a liquid handling (LiHa) arm with fixed tips, a robotic manipulator (RoMa) arm and a bespoke tilting rack with a Thermo Scientific Variomag® Teleshake. The tilting rack allows complete liquid uptake from large diameter wells, for example from a 6-well plate, as well as shaking for mixing and cell detachment during cell dissociation (Figure 2).



Figure 1 System overview; 1. Biosafety cabinet, 2. LiHa, 3. Te-Shake™ with chamber slides, 4. Reagent troughs and LiHa wash station, 5. Sterile system liquid, 6. RoMa, 7. Flask Flipper, 8. Te-Shake with tubes, 9. Tilting rack, 10. Centrifuge, 11. Liquid waste. Not all components are required for this application.



Figure 2 Tilting rack with shaker, 6-well and 24-well plates.

Cell lines

Four cell lines – HeLa, SaOs, CHO and HEK – with various adherence properties were selected (Table 1 and Figure 3).

Cell line	Definition	Adherence property
HeLa	Derived from human cervical cancer (Henrietta Lacks, 1951)	Medium adherence Tend to grow irregularly
SaOs	Sarcoma osteogenic cells	Strong adherence Regular growth
CHO	Chinese hamster ovary cells. Widely used for protein production	Low adherence Can grow in suspension
HEK	Human Embryonic Kidney cells. Transformed cells	Low adherence

Table 1 Cell lines used for experiments.




Figure 3 Microscopy images of used cell lines, x40.

Experimental conditions

Cells were grown in standard 6-well and/or 24-well plates (growth area 57 and 45.6 cm² respectively) until approximately 70 % confluency was reached. At this confluency cells are in the exponential growth phase and may be split or used for assays.

Optimal harvesting parameters, such as reagent volume, incubation time and mixing, were initially determined manually for each cell type. The protocol was subsequently transferred to the Freedom EVO liquid handling platform and the automated procedure optimized (Table 2).

Experiments were performed in parallel by both automatic and manual procedures, using three to six microplates. Cells from all the wells on a plate were combined and the overall cell number, viability and aggregation rate determined using a Cedex[®] (Roche Innovatis) cell counter. After harvesting, plates were checked under a microscope to confirm that all cells had been detached and collected.

	SaOs	Hela	Hela	CHO	HEK	HEK
Plate format	6 well	6 well	24 well	6 well	6 well	24 well
Growth medium with FCS 10%; P/S 1%	DMEM/F-12; 5 mM Glu	DMEM/F-12	DMEM/F-12	Ham's F-12, W 1.176 g/l, NaHCO ₃	DMEM/F-12	DMEM/F-12
Medium volume (ml)	2.0	2.0	0.8	3.0	2	0.8
PBS rinse (ml)	1.8	1.8	0.8	1.8	1.8	0.8
Volume Trypsin (ml)	1.2	1.2	0.5	1.0	1.0	0.4
Incubation (min)	I) 15, II) 15	15	16	5	3	5
Shaking (rpm)	1200	1200	1600	1200	1600	1200
Medium stop (ml)	1.2	1.2	0.6	1.2	1.6	0.6
Mix with LiHa (ml)	2.5	2.5	1.5	2.4	3	1.4
Medium rinse (ml)	1.5	1.5	0.6	1.2	1.5	0.6
Duration (min)	44	30	35	20	17	26


Table 2 Parameters for harvesting of cell lines in multiwell plates.

Results and Discussion


Automated harvesting was comparable to manual harvesting with regard to cell count, viability and aggregation rate (Figure 4).

Cell count for manual harvesting was standardized at 100 %, and results from the automated procedure measured relative to this value. Automated cell harvesting achieved cell counts ranging between 92-105 %, with viability between 80 % and 100 %. However, HeLa cells harvested from 24-well plates with a confluency of 90-100 % showed a viability of only 60 %. This low viability is probably the result of cells reaching the stationary stage due to nutrient depletion. Aggregation rate was highest for HEK cells, however the aggregation was not due to insufficient mixing as this cell type is known to build aggregates. Detailed results are shown in Figure 4.


SaOs (n = 4), 6-well plates




HeLa (n = 6), 6-well plates




HeLa (n = 4), 24-well plates, confluency 50-60 %




HeLa (n = 3), 24-well plates, confluency 90-100 %



CHO (n = 4), 6-well plates



HEK (n = 4), 6-well plates



HEK (n = 5), 24-well plates




Figure 4 Cell amount, cell viability and aggregation rate of manual and automated protocols in multiwell plates.

The most critical parameters for a successful cell harvest with sufficient density and viability were found to be confluency, dissociation procedure and pipetting speed. Protocols were optimized for a confluency of about 70 %, unless otherwise stated. An automated confluency check can assist determination of optimal harvesting time, and Tecan has integrated the Cellavista[®] cell imaging system – which allows for non-invasive brightfield and fluorescence analysis of cell culture flasks and microplates – onto the Freedom EVO platform.

Tapping the microplate to release dissociated cells is another important consideration when transferring a manual procedure to an automated system. Although the Teleshake was able to suspend cells, parameters such as shaking frequency and time needed to be optimized for each cell line. The speed was adjusted so that the entire well surface was covered with reagent and agitated during incubation. The dissociation reaction was stopped by the addition of medium containing 10 % serum and subsequent shaking step increased harvesting rate significantly (data not shown). The pipetting speed of the cell suspension, and uptake of partition volume, was very important to obtaining comparable harvesting results. Pipetting speed must be moderate, even if this prolongs the procedure for large volumes, as high pipetting speeds and extensive mix cycles might lead to a lower number of cells harvested. Good maintenance practice must be observed, with regular thorough cleaning of the liquid system even if it has not been used, in order to obtain reproducible results.

Conclusions

Comparable cell harvesting results were obtained for cell count, viability and aggregation rate for manual and automated procedures performed in multiwell plates. Individual protocols have been developed and automation parameters optimized for each of the four cell lines selected; this process must be repeated and adapted for additional cell lines. The Freedom EVOware scripts developed for this study, which are available from Tecan, provide a good starting point for protocols for new cell lines. The four cell lines studied were chosen for their wide range of characteristics, such as strength of adherence and aggregation rate, so the development of protocols for other cell lines with similar characteristics will be straightforward.

Acknowledgements

We would like to thank Prof Dr Ursula Graf-Hausner and Dr Nicola Francini, Institute of Chemistry and Biological Chemistry, Zurich University of Applied Sciences, Wädenswil, Switzerland, for the collaboration and for performing the experiments and establishing the protocols.

Austria +43 62 46 89 33 **Belgium** +32 15 42 13 19 **China** +86 21 289 8633 **Denmark** +45 70 23 44 50 **France** +33 4 72 76 04 80 **Germany** +49 79 51 94 170
Italy +39 02 92 44 790 **Japan** +81 44 556 73 11 **Netherlands** +31 18 34 48 174 **Singapore** +65 644 41 886 **Spain** +34 93 490 01 74
Sweden +46 31 75 44 000 **Switzerland** +41 44 922 89 22 **UK** +44 118 9300 300 **USA** +1 919 361 5200 **Other countries** +41 44 922 8125

Tecan Group Ltd. makes every effort to include accurate and up-to-date information within this publication, however, it is possible that omissions or errors might have occurred. Tecan Group Ltd. cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. All mentioned trademarks are protected by law. For technical details and detailed procedures of the specifications provided in this document please contact your Tecan representative. This brochure may contain reference to applications and products which are not available in all markets. Please check with your local sales representative.

Tecan, Freedom EVO and Freedom EVOware are registered trademarks and Te-Shake is a trademark of Tecan Group Ltd., Männedorf, Switzerland. RoboFlask is a registered trademark of Corning Incorporated, Corning, USA. Variomag is a registered trademark of H + P Labortechnik AG, Oberschleißheim, Germany. Cellavista and Cedex are registered trademarks of Roche Diagnostics GmbH, Mannheim, DE.

© 2010, Tecan Trading AG, Switzerland, all rights reserved.

www.tecan.com

